

Nepenthe Laboratory 12 panel Enzyme Immunoassay

For providers who do not perform presumptive screening in their office, Nepenthe offers a 12 panel EIA screen with confirmation of positive results, as well as unexpected negatives of submitted prescribed medications.

Drug/Drug Class	Lab EIA Cut Off	LC-MS/MS Cut Off⁵
Alcohol (Ethyl Glucuronide)	500ng/mL	250-500ng/mL
Amphetamines	500ng/mL	50-100ng/mL
Barbiturates	200ng/mL	100ng/mL
Benzodiazepines	200ng/mL	50ng/mL
Buprenorphine	20ng/mL	20-50ng/mL
THC	50ng/mL	50ng/mL
Cocaine	150ng/mL	50ng/mL
MDMA	500ng/mL	50ng/mL
Methadone	300ng/mL	50ng/mL
Opiates	300ng/mL	20-50ng/mL
Oxycodone	100ng/mL	50ng/mL
Tricyclic Antidepressants	300ng/mL	20-50ng/mL